

Case Report

DETECTION AND
VISUALIZATION OF
NEUROCENTRAL
SYNCHONDROSIS
WITH DSG CONNECT

Heiko Koller, PD Dr. med.

Facharzt für Orthopädie und Unfallchirurgie
at German Scoliosis Center, Werner-Wicker-
Clinic, Bad Wildungen, Germany

SpineGuard[®]
Making spine surgery safer



Case

8.5 years old patient (Fig.1) with progressive Early-Onset scoliosis;
Major thoracic curve
Risser 0-



Fig. 1

Treatment

Posterior transmuscular instrumentation with trolley-type growth guiding construct.

Pedicle Preparation: PediGuard® DSG Connect

The PediGuard DSG Connect probe is enabled with DSG® (Dynamic Surgical Guidance) Technology. It's a stand-alone, pedicle preparation instrument. A bipolar sensor is embedded at the tip of each PediGuard device. A low frequency and low voltage current is emitted through a bipolar electrode in the probe. Measuring the local conductivity of tissue 5 times per second, the device can accurately inform the surgeon of the type of tissue at its tip by changes in the pitch and cadence of an audio signal and a flashing LED light. Thus, the surgeon can be alerted of an imminent cortical breach and redirect appropriately during pedicle preparation.

Information about DSG Connect to add.



Fig. 2

Fig. 3

The PediGuard 2.5 XS probe was used during this case.

The PediGuard probe is used straight forward in manual techniques as the surgeon is used to with his standard pedicle awl.

Pedicle Preparation Technique

All of the pedicles were cannulated using a DSG Connect PediGuard 2.5 XS probe and free-hand technique.




Fig. 4



Fig. 5

Postoperative Follow-up

Fig. 6





DSG Technology:
Simple.
Accurate.
Smart.

SpineGuard®
Making spine surgery safer

www.spineguard.com

SpineGuard® S.A.
10, Cours Louis Lumière
94300 Vincennes, France
Phone: +33 1 45 18 45 19
Fax: +33 1 45 18 45 20

SpineGuard® Inc.
1434 Spruce Street, Suite 100
Boulder, CO 80302
Phone: +1 720 726 8085
Fax: +1 720 726 8001

CE 0123

Caution: See package insert for labeling limitations, intended uses, relevant warnings, precautions, side effects and contraindications. Federal (USA) law restricts the sale and use of this device to a prescription of a physician.