EU version:

EN - Instructions for use	2
FR - Notice d'utilisation	3
DE - Gebrauchsanweisung	4
ES - Instrucciones de uso	5
EL - Οδηγίες χρήσης	6
IT - Istruzioni d'uso	7
NL - Gebruiksaanwijzing	8
PT - Instruções de uso	9


U.S. version:

EN - Instructions for use	10
----------------------------------	----

FFX® IMPLANT (Facet FiXation)
INSTRUCTIONS FOR USE
DESCRIPTION

FFX® IMPLANTS are titanium devices (Ti-6Al-4V, ISO 5832-3) to be fixed in the facet joint in adult patients. They are implanted between the two facets to be stabilized and can contain bone graft or bone substitute. They can be fixed with or without transfacet screws. Properly impact implants between facets to avoid implant migration.

Reference	A (mm)	B (mm)	C (mm)	Color
57.000.10S	10	11	2.5	Gold
57.000.20S	10	13	2.5	Dark blue
57.000.30S	10	11	3	Green
57.000.40S	10	13	3	Light blue
57.000.50S	10	11	3.5	Brown
57.000.60S	10	13	3.5	Purple

INDICATIONS

FFX® IMPLANTS are intended for the treatment of the following conditions:

- Degenerative lumbar spinal stenosis
- Facet syndromes (primary or secondary)

INTENDED USE

FFX® implants prevent spinal instability and facet motion by enabling lumbar facet joint arthrodesis consecutive to the surgical treatment of degenerative lumbar spinal stenosis or facet syndrome.

CONTRA-INDICATIONS

The use of **FFX® implants** is contra-indicated in the following situations:

- Severe or chronic infections, local or systemic
- Sensitivity to materials
- Bone destruction or demineralization potentially affecting fixation of the implant
- Severe muscular, neurological, or vascular deficiencies affecting the limb in question
- Any infection that could compromise the function of the implant
- Morbid obesity of the patient
- Psychiatric background
- Osteoporosis (BMD < 120 mg/cm³)
- Major spine instability (spondylolisthesis ≥ grade II)
- Fracture of the spine, isthmus lysis
- Important scoliotic deformities (>25°)
- More than 3 vertebral levels concerned
- Spondylodiscitis or spine tumor
- Wide resection of facets joints during the surgery
- Unilateral application of the device

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GENERAL INSTRUCTIONS FOR USE
Pre-operative planning:

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- Improper positioning of an implant could result in its failure. Implants should therefore be placed only by a surgeon with the necessary training
- Do not use implants as test devices
- Dispose of implants and waste coming into contact with patients in accordance with current regulations
- Upon surgeon intraoperative assessment, if significant instability caused by decompression is observed (e.g., extended laminectomy), it is recommended to fix the FFX® implants with transfacet, auto-compressive screws made of titanium, of a length between 14 and 30 mm, and a diameter between 2 and 4.5 mm, depending on the patient's anatomy
- In the event of the sterile packaging being damaged or unintentionally opened before use, the device must be discarded
- The implant is a single use device, do not reuse

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POST-OPERATORY PRECAUTIONS

- Advise the patient of the precautions to be followed as part of post-operative development
- Carry out radiological surveillance (frequency and protocol determined by the surgeon)
- Proceed with timely and effective treatment of any infection (nearby cutaneous, urinary, bronchopulmonary, dental, etc.), however benign, due to risk of haematogenous contamination
- Be attentive to any signs of pain
- Consult our company if revision surgery appears necessary

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STERILITY

FFX® implants are delivered sterile. Sterilization is carried out by gamma radiation.

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MEANING OF SYMBOLS

	Catalogue number		Use-by date		Consult instructions for use
	Sterilized using irradiation		Lot number		Do not reuse
	Manufacturer		CE mark		Date of manufacture
	Do not use if package is damaged		Medical device		Double sterile barrier system

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INSTRUCTIONS FOR USE

DESCRIPTION

FFX® is a titanium device (Ti-6Al-4V, ISO 5832-3) to be fixed in the facet joint in adult patients. It is implanted between the two facets to be stabilized and must be packed with autograft. It must be fixed with transfacet screws. Properly impact the implant between facets to avoid implant migration

Reference	A (mm)	B (mm)	C (mm)	Color
57.000.10S	10	11	2.5	Gold
57.000.20S	10	13	2.5	Dark blue
57.000.30S	10	11	3	Green
57.000.40S	10	13	3	Light blue
57.000.50S	10	11	3.5	Brown
57.000.60S	10	13	3.5	Purple

INDICATIONS

FFX is a lumbar facet device that is placed bilaterally through a posterior surgical approach and spans the facet interspace. FFX must be used with an FDA-cleared transfacet screw cleared for use in the lumbar spine. FFX is intended to provide temporary fixation and stabilization to the spine as an aid to lumbar fusion through bilateral immobilization of the facet joints at one or two levels with autogenous and/or allogenic bone graft. FFX must be accompanied with an FDA-cleared intervertebral body fusion device and/or with an FDA-cleared posterior lumbar pedicle screw and rod system implanted at the same spinal level(s) as an adjunct to a single or two-level intervertebral body or posterolateral fusion, respectively. FFX is indicated for the treatment of patients with lumbar degenerative disc disease (DDD) from L3 to S1 in skeletally mature patients who have failed conservative care

CONTRA-INDICATIONS

The use of **FFX®** is contra-indicated in the following situations:

- Severe or chronic infections, local or systemic
- Allergy to titanium
- Bone destruction or demineralization potentially affecting fixation of the implant
- Severe muscular, neurological, or vascular deficiencies affecting the limb in question
- Any infection that could compromise the function of the implant
- Active or prior infection at the surgical site
- Morbid obesity of the patient (e.g., BMI > 40)
- Psychiatric background
- Osteoporosis (T score of -2.5 or lower)
- Major spine instability (spondylolisthesis ≥ grade II)
- Fracture of the spine, isthmus lysis
- Scoliotic deformities (Cobb Angle >25°)
- More than 3 vertebral levels requiring surgical treatment
- Non-contiguous lumbar levels requiring surgical treatment
- Spondylodiscitis or spine tumor
- Wide resection of facets joints during the surgery or other type of decompression requiring more bone removal than a laminectomy
- Unilateral application of the device

ADVERSE EFFECTS AND POSSIBLE COMPLICATIONS

The patient must be warned of the limitations and risks inherent to the placing of the implant. Some complications can require a revision surgery:

- Early or late-onset infections, haematoma, allergies, cutaneous necrosis
- Necrosis of bone or bone resorption, pseudarthrosis, chondrolysis, algodystrophy, bone resorption or fractures, consolidation delays
- Mechanical complications (implant failure or migration, wear, loosening, etc.), implant fracture
- Neurological complications (clinical or electrical monitoring is advised: wake-up test), pain
- Cardiovascular disorders, thrombosis

GENERAL INSTRUCTIONS FOR USE

Pre-operative planning:

The choice of the size of the implant to be used depends on patient pathology and morphology.

Surgical technique: Implants should be placed using the ancillary material provided for this purpose. Implants should be placed only by a surgeon with the necessary training. Videos and documents of the surgical techniques are available on SC MEDICA website, on demand.

Warnings:

- Do not use a stainless-steel screw with a titanium implant, to prevent any risk of corrosion
- Perform a visual inspection of the implant to detect any possible damage (scratches, marks, etc.). Do not allow the implant to come into contact with instruments that could mark its surface. Do not use an implant if it has been subject to any impact or has visible damages
- Our company cannot be held liable for any corrective modification made to the implants

- Do not use a stainless-steel screw with a titanium implant, to prevent any risk of corrosion
- Perform a visual inspection of the implant to detect any possible damage (scratches, marks, etc.). Do not allow the implant to come into contact with instruments that could mark its surface. Do not use an implant if it has been subject to any impact or has visible damages
- Our company cannot be held liable for any corrective modification made to the implants
- Improper positioning of an implant could result in its failure. Implants should therefore be placed only by a surgeon with the necessary training
- Do not use implants as test devices
- Dispose of implants and waste coming into contact with patients in accordance with current regulations
- FFX must be fixed with FDA-cleared transfacet, fully threaded screws made of titanium, of a length between 14 and 25 mm, and a diameter between 3.5 and 4 mm, depending on the patient's anatomy
- In the event of the sterile packaging being damaged or unintentionally opened before use, the device must be discarded
- The implant is a single use device, do not reuse

Some implants are made of multiple components. Do not combine components from different manufacturers or origins. In this event, our company waives all liability.

POST-OPERATORY PRECAUTIONS

- Advise the patient of the precautions to be followed as part of post-operative development
- Carry out radiological surveillance (frequency and protocol determined by the surgeon)
- Proceed with timely and effective treatment of any infection (nearby cutaneous, urinary, bronchopulmonary, dental, etc.), however benign, due to risk of hematogenous contamination
- Be attentive to any signs of pain
- Consult our company if revision surgery appears necessary

The surgeon is to decide on authorizing free movement of the patient and to define any limits on patient activity after implantation. However, excessive activity involving the affected limb is strongly discouraged.

INTERFERENCE WITH OTHER MEDICAL TESTS OR TREATMENTS

- Intramuscular injections on the same side as the affected limb are strongly discouraged
- The presence of titanium in the implant generates artefacts in the immediate vicinity of the device. This must be taken into account in relation to medical imaging tests
- In case of breakage of an instrument during implantation, both the instrument and the implant must be replaced

HANDLING-STORAGE

Follow the required precautions for asepsis instructions when removing from packaging. Some products have sharp components that could injure the handler. The devices being provided sterile, their packaging is not to be opened until time of use. Store products in their packaging, in a dry place, out of direct sunlight and at room temperature.

STERILITY

FFX® is delivered sterile. Sterilization is carried out by gamma radiation.

Re-sterilization of products by the user is not allowed.

Warning and precaution: the implant is a single-use device and must not be reused. Mechanical, physical or chemical characteristics of the device may be altered if it is repeatedly cleaned and sterilized, thereby compromising the integrity of the structure and / or material of the device and, hence, the safety, performance and / or compliance with current standards.

MR1 SAFETY INFORMATION

Non-clinical testing has demonstrated that the **FFX** range is MR Conditional in accordance with the ASTM F2503-23e1 standard definitions.

MR1 Safety Information

A person with the SC Medica FFX implant may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.

Device Name	FFX® (Facet FiXation)
Static Magnetic Field Strength (B ₀)	1.5T or 3.0T
Maximum Spatial Field Gradient	30 T/m (3,000 gauss/cm)
RF Excitation	Circularly Polarized
RF Transmit Coil Type	Whole-Body Transmit RF Coil
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	2W/kg (Normal Operating Mode)
Scan Duration	2 W/kg whole-body average SAR for 60 minutes of continuous RF (a sequence or back to back series/scan without breaks)
MR Image Artifact	The presence of this implant may produce an image artifact.

ABLATION AND DISPOSAL OF DEVICE

The surgeon makes the final decision concerning ablation of an implant. Disposal must be carried out according to the procedures in force in the hospital establishment.

MEANING OF SYMBOLS

	Catalogue number		Use-by date		Consult instructions for use
	Sterilized using irradiation		Lot number		Do not reuse
	Manufacturer		Rx only		Date of manufacture
	Do not use if package is damaged		Medical device		Double sterile barrier system
	MR conditional				

For additional information, please contact the company SC MEDICA.

Any serious incident that has occurred in relation to the device should be reported to SC MEDICA and the competent authority of your State.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician. For Professional Use Only.

SC MEDICA
3 Quai Kléber, Tour Sébastopol
67000 Strasbourg – France
www.sc-medica.com
Tel. : +33 3 88 23 71 04
Fax: +33 3 67 10 97 53

Document # : DEV-REC-User manual_FFX_V11/ 09.2024

FFX[®]

Facet FiXation