FR - Notice d'utilisation	2
EN - Instructions for use	3
DE - Gebrauchsanweisung	4
ES - Instrucciones de uso	5
EL - Οδηγίες χρήσης	6
IT - Istruzioni d'uso	7
NL - Gebruiksaanwijzing	8
PT - Instruções de uso	9

FFX[®]

Instruments

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INSTRUCTIONS FOR USE

DESCRIPTION AND INTENDED PURPOSE

The **FFX® implant ancillary instruments** are a set of reusable surgical instruments, exclusively dedicated to the implantation and the explantation (if applicable) of **FFX® IMPLANTS**. The ancillary consists of the instrumentation and a container.

Description	Material	Expected performance
Facet-holder	Stainless steel	Hold FFX® and position it.
Targeting devices / working tubes	Titanium Stainless steel	Aligning the facet screw with the middle of FFX® implant
Impactor	Stainless steel	Drive the FFX® implants until there is enough posterior behind to graft.
Facet chisel	Stainless steel Titanium	Track the articular line spacing
Periosteal Elevator	Stainless steel	Track the curved articular line spacing
Rasp	Stainless steel Titanium	Expose the articular facets before placing the FFX implants
Graft holder	Titanium	Implant is placed into the appropriate place in the graft holder to fill the empty space with the filling material (bone graft or bone substitute)
Implant trial FFX 11x2.5/3	Titanium	Test the space at facet level and select the appropriate implant size.
Interlaminar distractor	Stainless steel	Open the facet space
Curette	Stainless steel	Bone scraping
Mallet	Stainless steel, Titane, Propylux	Insert/extract the instruments
Screwdriver shaft	Stainless steel	Insert the screw in the facet bone
Drill	Stainless steel	Drill the facet bone before screw implantation
Awl	Stainless steel	Make a pre-hole in the bone

INITIAL USE AND REPROCESSING

The ancillaries are delivered non-sterile, with a dedicated container. This section aims to provide the necessary instructions for the initial treatment and subsequent reprocessing of the ancillaries. All ancillary equipment must imperatively be decontaminated, cleaned and sterilized before use by the healthcare establishment, according to validated methods. The recommendations below do not replace the health rules in force (good practices, national recommendations, regulations, etc.). SC MEDICA performs routine inspections of ancillary equipment upon receipt and distribution. The lifetime of non-sterile instruments depends on several factors including, but not limited to, the method and duration of use and handling of the instruments. SC MEDICA recommends careful functional inspection of each instrument, to ensure that no signs of wear / damage are present and to ensure their effectiveness.

Device inspection

Upon receipt, inspect the packaging and labeling of the product. Upon receipt, and before each use, check the integrity of the product, in particular the absence of corrosion, scratches, nicks, residue, debris or change in color. after unpacking. In case of damage, the product must be disposed of. The mallet is fitted with a plastic tip which must be replaced as soon as signs of wear appear

Automatic cleaning:

Initial / point-of-use treatment

Before first use and all subsequent reuses, instruments made up of dismantlable components must be disassembled before pre-disinfection and cleaning: unscrew the implant holder, unscrew the plastic nozzle of the mallet, disassemble the screwdriver.

After each use, remove soils/stains with absorbent paper or a disposable lint-free wipe before they dry. Flush cannulated instruments several times with deionized water or distilled water.

Caution: the curette has a sharp blade. Handle it with care and treat it with a brush to avoid any risk of injury.

Immerse the instruments in an enzymatic solution of alkaline pH or neutral pH, specially designed for surgical instruments (for example, ANIOSYME X3 - 0.5%), for minimum 15 min. at room temperature. Follow the recommendations of the solution's manufacturer for dilution, temperature, and soak time. Make sure that all surfaces of the instrument are in contact with the solution and that no air remains trapped in any corners or holes.

Use a soft bristle brush to clean the instrument, paying particular attention to rough areas, corners and holes. Rinse the instrument by agitation under the water level (raw water) for at least two minutes.

Perform a visual inspection to make sure that all soiling has been removed, failing which, repeat the previous steps.

Washer-disinfector

SC MEDICA recommends an automatic cleaning, using a washer-disinfector, compliant with EN ISO 15883-1 standard :

Steps	Product	Duration	Temp.	Water
Precleaning	Water	2 min.	Room temp.	Raw water

ENGLISH

Cleaning	Alkaline detergent with surfactant	5 min.	55°C	Raw water
Rinsing	Water	2 min.	Room temp.	Purified Water
Thermal disinfection	Rinsing liquid for automated processing	6 min.	93°C	Purified Water
Drying	/	30 min.	Air : 110°C (min)	/
Cooling	/	5 min.	Air : 30°C	/

An alkaline or neutral cleaning product (pH: 6.0-8.5, e.g., NEODISHER MEDICLEAN 0.2-1%) must be used.

Follow the recommendations of the manufacturer of the washing and disinfecting solutions regarding dilution, temperature and duration of soaking. In case of insufficient drying, dry the implant using a blower (medical air).

Inspection

A visual inspection should be performed prior to sterilization. All instrument parts should be inspected to verify that all traces of soiling have disappeared, paying particular attention to the corners of the instruments. Do not use a damaged instrument. Proper drying of instruments prior to packaging for sterilization and storage is essential, as a wet instrument may compromise sterilization.

Sterilization

Place the instruments in the sterilization tray. Orient the hollow part of the mallet head towards the bottom of the basket.

SC MEDICA recommends performing the following sterilization process :

Temperature	134°C
Exposure time	18 minutes
Drying time	20 minutes
Method	Moist heat sterilization according to ISO 17665-1
Cycle	Saturated steam with forced fractional air removal

Follow the loading preparation instructions provided by the sterilizer manufacturer.

Storage

After sterilization, the instruments should be stored in the sterilization packaging in a dry and dust-free place.

Return after use

After each use and before returning the devices to SC MEDICA, the ancillary must be pre-disinfected, cleaned and sterilized according to the above recommendations. Any non-functional instrument must immediately be returned to SC MEDICA for maintenance or exchange. The nature of the malfunction must be reported.

The instrumentation must be returned properly packed in the corresponding tray. It must be accompanied by its traceability sheet duly completed by an authorized signatory by the healthcare facility.

MATERIAL DISPOSAL

Equipment disposal must be carried out according to the procedures in force in the healthcare facility.

WARNINGS

- The surgeon must be familiar with FFX® implants, the method of application, the instruments and the operating technique.
- Under no circumstances must ancillaries for FFX® implants be combined with other ranges of devices.
- It is recommended to check the allergic level of patients with particularly strong ground. Do not expose patients allergic to titanium and steel (nickel) to FFX ancillaries.
- Before each use, carry out a careful visual examination of the instrument in order to detect any possible deterioration (scratches, pitting, marking defects, etc.). In the event of a visible defect or damaged material, do not use the instrument. Microscopic defects due to reuse can cause rupture during surgery.
- Do not put the devices in contact with instruments that could alter their surface.
- As soon as an anomaly is noted, contact the manufacturer in order to take the appropriate measures.
- Verify the compatibility of the navigation instruments with the navigation system of the health facility.

MEANING OF SYMBOLS

	Catalogue number		Non-sterile		Consult instructions for use
	CE mark		Lot number		Manufacturer
	Medical device		Date of manufacture		Unique Device Identifier

Year of first CE marking: 2021

For additional information, please contact the company SC MEDICA. Any serious incident that has occurred in relation to the device should be reported to SC MEDICA and the competent authority of your State.

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