



OGMEND®

· Implant
· Enhancement
· System



**STRONG, FAST, SIMPLE
SOLUTION TO ENHANCE FIXATION
IN COMPROMISED SCENARIOS.**



A variety of scenarios make it difficult to achieve or maintain high-quality fixation between hardware and bone during instrumented spine surgery.



THE OGMEND® ADVANTAGE

Insert an Ogmend® Sleeve prior to inserting a screw to achieve stable fixation by enhancing the screw-to-bone interface



Implant in less than 2 minutes




Cost efficient for immediate savings



Used with preferred fixation systems



Sterile packed with no logistics needed



Accommodates wide range of screw sizes




Positive Data


MECHANISMS OF ACTION






Trap screw for immediate increased fixation strength.




Reduce and distribute stress to enhance the interface and reduce overload.



Enable strong, rapid bony ingrowth to maintain effect over time (does not resorb).





USE CASES

Ogmend® applies to a broad range of surgeries in the following scenarios:

Pre-Op Anticipated


Predicted compromised purchase based on pre-operative assessment

Intraoperatively

Insufficient purchase or loss of fixation during index procedure


Revision/Extension

During secondary or revision surgeries to correct failed screws



Simply...

Cut to Size



Insert into Bone



Insert Screw




Physicians implanting the device must refer to the instructions for use for a complete listing of procedures, contraindications, warnings, and precautions associated with the use of the Ogmend® Implant in their respective region.



ACCESS

Ogmend® maintains
FDA and CE Mark clearance for
use in spine fixation surgery.



Currently Available Sizes:

- Medium: Accommodates screw diameters from 3.5mm - 6.5mm
- Large: Accommodates screw diameters from 6.5mm - 10.5mm*

To Place an Order

- Contact sales: sales@wovenorthopedics.com
- Order online: wovenortho.com/fda or wovenortho.com/cemark



*Large is currently available only to countries who accept CE Mark.