

pediatric plates

proximal femur





Indications

- Developmental Dysplasia of the Hip, Cerebral Palsy,
- Legg-Calve Perthes Disease, Slipped Capital Femoral Epiphysis.

Surgical applications in this context can be listed as:

- Intertrochanteric osteotomy,
- Varus osteotomies with excision of the trochanter minor,
- Valgus and valgus rotation osteotomies for the femoral heads,
- Rotation osteotomies above or below the trochanter minor as required.





Varus 100°

Varus
110°

Varus
120°

Valgus
150°

Trauma Plate
120°



2-5 years of age and up to 20 kg → 2.7 mm Small Plates



2-9 years of age and up to 40 kg → 3.5 mm Small Plates

pediatric plates


distal femur

Indication

- Distal femur derotation and varus osteotomies.

Features

- Usable for both left and right femur.
- Angular stability reduces the risk of primary and secondary loss of correction.



pediatric plates

distal tibia



Indication

- Pediatric tibia fractures and osteotomies.

Features

- Usable for both left and right tibia.
- Designed with an angled-extended metaphyseal part
- Low profile in order to reduce skin irritation.





9-17 years of age and up to 60 kg → 4.5 mm Standard Plates

